
AIM-CU

Release 1.0.0

Smriti Prathapan, Berkman Sahiner, Dhaval Kadia, Ravi K. Samala

Jul 25, 2025

CONTENTS:

1	System setup	2
2	Code execution	3
3	Example code execution	4
4	Demo	5
5	Usability	6
6	AIM-CU	7
6.1	Methods	7
6.2	CUSUM	8
6.3	ARLTheoretical	9
6.4	Utils	10
7	Disclaimer	12
	Python Module Index	13
	Index	14

A CUSUM-based tool for AI Monitoring

Monitoring a clinically deployed AI device to detect performance drift is an essential step to ensure the safety and effectiveness of AI.

AIM-CU is a statistical tool for AI monitoring using cumulative sum (AIM-CU).

AIM-CU computes:

- The parameter choices for change-point detection based on an acceptable false alarm rate
- Detection delay estimates for a given displacement of the performance metric from the target for those parameter choices.

SYSTEM SETUP

Make sure R is installed in the system. Instructions for linux (the below setup is only performed in linux):

```
wget -q0- https://cloud.r-project.org/bin/linux/ubuntu/marutter_pubkey.asc | tee -a /  
↳etc/apt/trusted.gpg.d/cran_ubuntu_key.asc  
add-apt-repository "deb https://cloud.r-project.org/bin/linux/ubuntu $(lsb_release -cs)-  
↳cran40/"  
apt-get install -y --no-install-recommends r-base r-base-dev  
  
# setup R configs  
echo "r <- getOption('repos'); r['CRAN'] <- 'http://cran.us.r-project.org';  
↳options(repos = r);" > ~/.Rprofile  
Rscript -e "install.packages('ggplot2')"  
Rscript -e "install.packages('hexbin')"  
Rscript -e "install.packages('lazyeval')"  
Rscript -e "install.packages('cusumcharter')"  
Rscript -e "install.packages('RcppCNPY')"  
Rscript -e "install.packages('spc')"
```

CODE EXECUTION

Clone AIM-CU repository.

```
git clone https://github.com/DIDSR/AIM-CU.git
```

Run the following commands to install required dependencies (Python = 3.10 is used).

```
apt-get -y install python3  
apt-get -y install pip  
cd AIM-CU  
pip install -r requirements.txt
```

Run AIM-CU.

```
cd src/package  
python3 app.py
```

Open the URL <http://0.0.0.0:7860> that is running the AIM-CU locally.

EXAMPLE CODE EXECUTION

Example code can be run in a Jupyter Notebook after opening it with `jupyter notebook` command from `/src/package/` directory. The tool is designed to be used through UI, not from console.

CHAPTER FOUR

DEMO

AIM-CU can also be run through the demo available at <https://huggingface.co/spaces/didsr/AIM-CU>. If Space is paused, click on Restart button.

USABILITY

- Example AI output CSV file is available as [config/spec-60-60.csv](#) to be uploaded in monitoring phase.
- Workflow instruction to run the tool is available at bottom-left of UI.
- Sample UI output is available at [assets/ui.png](#).
- Setting `control:save_figure` to `true` from [config.toml](#) will save tables and plots in [figure/](#).
- Running AIM-CU does not take time longer than a few seconds, and it does not require GPU.

6.1 Methods

6.1.1 CUSUM parameters

Table 1: CUSUM parameters

Pa- rame- ter	Description
μ_{in}	The mean of the performance metric when the process is in-control, i.e., when there is no performance drift
ARL_0	Number of observations before the control chart signals a false detection
σ_{in}	The in-control standard deviation of the metric
ARL_1	Number of observations before the control chart signals a true detection
k	The normalized reference value, which is related to the magnitude of change that one is interested in detecting. $k = 0.5$ is the default choice for detecting a unit standard deviation change
S_hi	Cumulative sum of positive changes in the metric
h	The normalized threshold or control limit (default =4). This threshold determines when the control chart signals a detection
S_lo	Cumulative sum of negative changes in the metric

6.1.2 CUSUM chart

A two-sided CUSUM control chart computes the cumulative differences or deviations of individual observations from the target mean (or in-control mean, μ_{in}). The positive and negative cumulative sums are calculated:

$$S_{hi}(d) = \max(0, S_{hi}(d-1) + x_d - \hat{\mu}_{in} - K)$$

$$S_{lo}(d) = \max(0, S_{lo}(d-1) - x_d + \hat{\mu}_{in} - K)$$

where d denotes a unit of time, x_d is the value of quantity being monitored at time d , $\hat{\mu}_{in}$ is the in-control mean of x_d , and K is a “reference value” related to the magnitude of change that one is interested in detecting. S_{hi} and S_{lo} are the cumulative sum of positive and negative changes. To detect a change in the observed values from the in-control mean, the CUSUM scheme accumulates deviations that are K units away from the in-control mean. Let σ_{in} denote the in-control standard deviation of x_d .

6.2 CUSUM

Cumulative Sum (CUSUM)

@author: smriti.prathapan

class package.cusum.CUSUM

CUSUM class and its functionalities.

change_detection(*normalized_ref_value: float = 0.5, normalized_threshold: float = 4*) → None

Detects a change in the process.

Parameters

- **pre_change_days** (*int*) – Number of days for in-control phase.
- **normalized_ref_value** (*float, optional*) – Normalized reference value for detecting a unit standard deviation change in mean of the process. Defaults to 0.5.
- **normalized_threshold** (*float, optional*) – Normalized threshold. Defaults to 4.

compute_cusum(*x: list[float], mu_0: float, k: float*) → tuple[list[float], list[float], list[float]]

Compute CUSUM for the observations in x

Parameters

- **x** (*list[float]*) – Performance metric to be monitored
- **mu_0** (*float*) – In-control mean of the observations/performance metric
- **k** (*float*) – Reference value related to the magnitude of change that one is interested in detecting

Returns

Positive cumulative sum, negative cumulative sum, and CUSUM

Return type

tuple[list[float], list[float], list[float]]

initialize() → None

Initialize with the configuration file.

plot_cusum_plotly() → Figure

Plot CUSUM value using Plotly

Returns

CUSUM plot using Plotly graph object.

Return type

go.Figure

plot_input_metric_plotly() → Figure

Plot the input metric using Plotly.

Returns

Scatter plot as Plotly graph object.

Return type

go.Figure

plot_input_metric_plotly_raw() → Figure

Plot AI output using Plotly.

Returns

Scatter plot as Plotly graph object.

Return type

go.Figure

set_df_metric_csv(*data_csv: DataFrame*) → None

Assign the performance metric data to be used for CUSUM.

Parameters

data_csv (*DataFrame or TextFileReader*) – A comma-separated values (csv) file is returned as two-dimensional data structure with labeled axes.

set_df_metric_default() → None

Read the provided performance metric data to be used for CUSUM for an example.

set_init_stats(*init_days: int*) → None

Use number of baseline observations to calculate in-control mean and standard deviation.

Parameters

init_days (*int, optional*) – Number of baseline observations when observations are considered stable. Defaults to 30.

set_timeline(*data: ndarray*) → None

Set the timeline of observations.

Parameters

data (*np.ndarray*) – Data of the metric values across the observations.

6.3 ARLTheoretical

ARLTheoretical

@author: smriti.prathapan

package.ARLTheoretical.get_ARL_1(*h: float, shift_in_mean: list[float], dict_ARL0_k: OrderedDict*) → DataFrame

Get the ARL1 along with k values.

Parameters

- **h** (*float*) – Normalized threshold.
- **shift_in_mean** (*list[float]*) – List of the values of shift in mean.
- **dict_ARL0_k** (*OrderedDict*) – Data dictionary of ARL0 and k

Returns

Table for ARL1 and k values.

Return type

pd.DataFrame

package.ARLTheoretical.get_ARL_1_h_mu1_k(*h: float, k: float, mu1: float*) → float

Calculate ARL_1 with given Shift in Mean (mu1) and k.

Parameters

- **h** (*float*) – Normalized threshold.
- **k** (*float*) – Normalized reference value.

- **mul** (*float*) – Intended shift in mean.

Returns

Detection delay (ARL1).

Return type

float

`package.ARLTheoretical.get_ref_value(h: float, list_ARL_0: list[float]) → tuple[DataFrame, OrderedDict]`
 provides normalized reference values k for provided list of ARL0, given the value of normalized threshold h.

Parameters

- **h** (*float*) – Normalized threshold.
- **list_ARL_0** (*list*) – List of ARL0 values.

Returns

Dataframe of ARL0 and k, Data dictionary of ARL0 and k; where k is normalized reference value.

Return type

tuple[pd.DataFrame, OrderedDict]

`package.ARLTheoretical.get_ref_value_k(h: float, ARL_0: float) → float`
 Calculation for the reference value for given h and ARL_0.

Parameters

- **h** (*float*) – Normalized threshold.
- **ARL_0** (*float*) – ARL0 value.

Returns

Normalized reference value k.

Return type

float

6.4 Utils

Utilities to handle different operations

`package.utils.get_greatable_as_html(df: DataFrame) → GT`
 Get the great_table as HTML from Pandas dataframe.

Parameters

df (*pd.DataFrame*) – Dataframe to render as a table.

Returns

Table in HTML format.

Return type

gt.GT

`package.utils.populate_summary_table_ARL0_k(summary_table_df_ARL0_k: DataFrame, h) → GT`
 Populate ARLTheoretical.summary_table_df_ARL0_k.

Parameters

- **summary_table_df_ARL0_k** (*pd.DataFrame*) – Dataframe of ARL0 and its respective values of k.

- **h** (*float*) – Normalized threshold.

Returns

Table of ARL0 and k in HTML format.

Return type

gt.GT

`package.utils.populate_summary_table_ARL1_k(summary_table_df_ARL1_k: DataFrame, dict_ARL0_k: OrderedDict, h) → GT`

Populate Multiindex table specific for ARLTheoretical.summary_table_df_ARL1_k

Parameters

- **summary_table_df_ARL1_k** (*pd.DataFrame*) – Dataframe with ARL1 and k values.
- **dict_ARL0_k** (*OrderedDict*) – Data Dictionary with the mapping between ARL0 and k.
- **h** (*float*) – Normalized threshold.

Returns

Table for ARL1 and k in HTML format.

Return type

gt.GT

DISCLAIMER

This software and documentation was developed at the Food and Drug Administration (FDA) by employees of the Federal Government in the course of their official duties. Pursuant to Title 17, Section 105 of the United States Code, this work is not subject to copyright protection and is in the public domain. Permission is hereby granted, free of charge, to any person obtaining a copy of the Software, to deal in the Software without restriction, including without limitation the rights to use, copy, modify, merge, publish, distribute, sublicense, or sell copies of the Software or derivatives, and to permit persons to whom the Software is furnished to do so. FDA assumes no responsibility whatsoever for use by other parties of the Software, its source code, documentation or compiled executables, and makes no guarantees, expressed or implied, about its quality, reliability, or any other characteristic. Further, use of this code in no way implies endorsement by the FDA or confers any advantage in regulatory decisions. Although this software can be redistributed and/or modified freely, we ask that any derivative works bear some notice that they are derived from it, and any modified versions bear some notice that they have been modified.

PYTHON MODULE INDEX

p

`package.ARLTheoretical`, [9](#)
`package.cusum`, [8](#)
`package.utils`, [10](#)

INDEX

C

`change_detection()` (*package.cusum.CUSUM method*), 8
`compute_cusum()` (*package.cusum.CUSUM method*), 8
`CUSUM` (*class in package.cusum*), 8

G

`get_ARL_1()` (*in module package.ARLTheoretical*), 9
`get_ARL_1_h_mu1_k()` (*in module package.ARLTheoretical*), 9
`get_greatable_as_html()` (*in module package.utils*), 10
`get_ref_value()` (*in module package.ARLTheoretical*), 10
`get_ref_value_k()` (*in module package.ARLTheoretical*), 10

I

`initialize()` (*package.cusum.CUSUM method*), 8

M

`module`
 package.ARLTheoretical, 9
 package.cusum, 8
 package.utils, 10

P

`package.ARLTheoretical`
 module, 9
`package.cusum`
 module, 8
`package.utils`
 module, 10
`plot_cusum_plotly()` (*package.cusum.CUSUM method*), 8
`plot_input_metric_plotly()` (*package.cusum.CUSUM method*), 8
`plot_input_metric_plotly_raw()` (*package.cusum.CUSUM method*), 8
`populate_summary_table_ARL0_k()` (*in module package.utils*), 10

`populate_summary_table_ARL1_k()` (*in module package.utils*), 11

S

`set_df_metric_csv()` (*package.cusum.CUSUM method*), 9
`set_df_metric_default()` (*package.cusum.CUSUM method*), 9
`set_init_stats()` (*package.cusum.CUSUM method*), 9
`set_timeline()` (*package.cusum.CUSUM method*), 9